### Application Bulletin 416/1 e

### Determination of arsenic in water with the scTRACE Gold

### Summary

Arsenic is ubiquitous in the crust of the earth in small concentrations. But mineral and ore deposits with elevated arsenic content can occur under special geological conditions (volcanic rock, phosphate or sulfide containing mineral deposits). From these deposits, arsenic is eluted as arsenite (AsO<sub>3</sub><sup>3-</sup>) and arsenate (AsO<sub>4</sub><sup>3-</sup>) and contaminates the ground water. In addition to this natural origin, also industry and agriculture can contribute to the contamination of water with arsenic as anthropogenic sources. The WHO (World Health Organization) recommends a maximum content of  $\beta$ (As) = 10 µg/L for water used as drinking water.

This Application Bulletin describes the determination of arsenic in water samples by anodic stripping voltammetry using the scTRACE Gold sensor. Furthermore, it is possible to distinguish between As(total) and As(III). With a deposition time of 60 s, the limit of detection for As(total) is  $0.9 \mu g/L$ , for As(III) it is  $0.3 \mu g/L$ .

### Samples

Surface water, ground water

### Instruments

797 VA Computrace	2.797.0020
Accessories	
Stirrer	6.1204.200
Driving belt	6.1244.020
Measuring vessel 5 mL	6.1415.150
SGJ Stopper B-14	6.1446.000
Сар	6.2753.210

### Electrodes

scTRACE Gold	6.1258.000
Electrode shaft	6.1241.080

**Note!** It is recommended to have a separate measuring vessel and stirrer for this application which have not been used in combination with a mercury or platinum electrode.

## The scTRACE Gold

The scTRACE Gold sensor holds all three electrodes required for a voltammetric determination. Together with the electrode shaft it makes a complete electrode system which can be used in any Metrohm voltammetric measuring stand. No further electrodes are required. The working electrode is a gold microwire. Reference and auxiliary electrode are screen printed electrodes.

Different from other solid state electrodes the scTRACE Gold does not need extensive conditioning before it can be used. A new sensor only needs to be activated as described in paragraph «Activation of the scTRACE Gold». This takes about 10 minutes and the sensor is ready for the first determination.

Due to its construction the scTRACE Gold is maintenancefree. It can be electrochemically cleaned as described in paragraph «Cleaning of the scTRACE Gold», but mechanical cleaning is neither necessary nor possible. As any electrode the performance of the scTRACE Gold will deteriorate with the number of determinations. Signals will get smaller and curves are less reproducible. Then it is time to replace the sensor.



### Reagents

- As(V) standard stock solution, β(As(V)) = 1 g/L, commercially available
- As(III) standard stock solution, β(As(III)) = 1 g/L, commercially available
- Sulfuric acid, w(H<sub>2</sub>SO<sub>4</sub>) = 96%, for trace analysis\*, CAS 7664-93-9
- Nitric acid, w(HNO<sub>3</sub>) = 65%, for trace analysis\*, CAS 7697-37-2
- Sulfamic acid,  $NH_2SO_3H$ ,  $\geq$  99.5%, analytical grade, CAS 5329-14-6
- Citric acid monohydrate, C<sub>6</sub>H<sub>8</sub>O<sub>7</sub> · H<sub>2</sub>O, for trace analysis\*, CAS 5949-29-1
- Potassium chloride, KCl, for trace analysis\*, CAS 7447-40-7
- Ultrapure water, resistivity >18 MΩ·cm (25 °C), type I grade (ASTM D1193)
- $^{\ast}$  e.g., Merck suprapur  $^{\otimes},$  Sigma-Aldrich TraceSelect  $^{\otimes}$  or equivalent

### Solutions

Cleaning solution	$c(H_2SO_4) = 0.5 \text{ mol/L}$ $c(KCI) = 0.05 \text{ mol/L}$ In a 100 mL volumetric flask 0.373 g KCI are dissolved in approx. 80 mL ultrapure water. 2.78 mL $w(H_2SO_4) = 96\% \text{ are carefully}$ added. Attention! Solution gets very hot. After cooling down to room temperature the solution is made up to the mark with ultrapure water.
Electrolyte	c(sulfamic acid) = 1 mol/L c(citric acid) = 0.5 mol/L c(KCI) = 0.45 mol/L 9.71 g sulfamic acid, 10.51 g citric acid and 3.35 g KCI are dissolved in 100 mL ultrapure water. An ultrasonic bath is used to dissolve everything.

Standard solutions	
As(V) standard	$\begin{split} \beta(As(V)) &= 1 \text{ mg/L} \\ Approx. 40 \text{ mL ultrapure water are} \\ filled into a 50 \text{ mL volumetric flask.} \\ 0.05 \text{ mL } w(HNO_3) &= 65\% \text{ and } 0.05 \\ \text{mL } As(V) \text{ standard stock solution} \\ \text{are added. The solution is made} \\ \text{up to the mark with ultrapure} \\ water. \end{split}$
As(III) standard	$\beta$ (As(III)) = 1 mg/L Approx. 40 mL oxygen-free ultrapure water are filled into a 50 mL volumetric flask. 0.05 mL w(HNO <sub>3</sub> ) = 65% and 0.05 mL As(III) standard stock solution are added. The solution is made up to the mark with oxygen-free ultrapure water. When the standard is stored dark and cool, it is stable for approx. 1 week.

# Activation of the scTRACE Gold

A new sensor needs to be activated. The activation only has to be carried out prior to the first use. If the electrode needs to be cleaned before, in between or after determinations the procedure described in the chapter «Cleaning of the scTRACE Gold» should be used.

### Analysis

10 mL of the cleaning solution are pipetted into the measuring vessel. The activation is carried out using the parameters given under «Parameters for activation».

### Measuring solution

10 mL cleaning solution

#### Parameters for activation

Determination	
No. of additions	0
No. of replications	4
Voltammetric	
Electrode	RDE/SSE

Application Bulletin 416/1 e	
Determination of arsenic with the scTRACE G	old

Measuring mode	DC – Sampled direct current
Stirring speed	2400 min <sup>-1</sup>
Hydrodynamic measurement	No
Conditioning cycles	
Start potential	-1.5 V
End potential	1 V
No. of cycles	10
Pretreatment	
Cleaning potential	0.1 V
Cleaning time	5 s
Deposition potential	-0.3 V
Deposition time	10 s
Sweep	
Equilibration time	5 s
Start potential	-0.3 V
End potential	0.2 V
Potential step	0.6 V
Potential step time	0.025 s
Sweep rate	24 V/s

### Example for the curves obtained during activation



### Comments

 The lifetime of the sensor will significantly be reduced if the activation is repeated or carried out more than once.

# Method 1: Determination of total arsenic

With method 1 the sum of As(III) and As(V) in the sample is determined. The As(V) species is electrochemically inactive. During the determination it is reduced in-situ by nascent hydrogen to As(III). Together with the As(III) present in the sample it is further reduced electrochemically to As(0) and deposited on the gold working electrode in the same step. During the subsequent stripping step the deposited As(0) is reoxidized to As(III) giving the analytical signal.

With 60 s for deposition the limit of detection is  $0.9 \ \mu g/L$  and the calibration is linear up to a concentration of 20  $\ \mu g/L$ .

### Sample preparation

Ground, drinking, and mineral water can usually be analyzed directly.

Water with a low to medium degree of contamination by organic substances is digested in the 909 UV Digester: 10 mL acidified water sample (pH = 2) are mixed with 10  $\mu$ L w(HNO<sub>3</sub>) = 65% and 50  $\mu$ L w(H<sub>2</sub>O<sub>2</sub>) = 30% and irradiated for 90 min at 90 °C.

### Analysis

10 mL sample are pipetted into the measuring vessel and 2 mL electrolyte are added. The determination is carried out using the parameters given under «Parameters for As(total)».

The concentration of As(total) is quantified by two additions of As(V) standard solution.

### Measuring solution

10 mL sample 2 mL electrolyte

#### Parameters for As(total)

Determination	
No. of additions	2
No. of replications	2
Voltammetric	
Electrode	RDE/SSE
Measuring mode	SqW – Square wave
Stirring speed	2400 min <sup>-1</sup>
Hydrodynamic measurement	No
Conditioning cycles	
Start potential	-0.2 V

End potential	1 V
No. of cycles	10
Pretreatment	
Cleaning potential	-1.0 V
Cleaning time	60 s
Deposition potential	(-0.25 V) see «Comments»
Deposition time	5 s
Sweep	
Equilibration time	5 s
Start potential	-0.3 V
End potential	0.4 V
Potential step	0.01 V
Amplitude	0.02 V
Frequency	100 Hz
Sweep rate	1 V/s
Substance + calibration	
Calibration	Standard addition
Name	As(total)
Peak potential	0.0 V
Tolerance	0.05 V
Baseline	Linear
	Automatic

### Example for the determination of As(total)



Application Bulletin 416/1 e Determination of arsenic with the scTRACE Gold



Results

Sample	β(As(total))
Bottled mineral water	4.9 μg/L

### Limit of detection and linear working range

The limit of detection was determined using the «regression approach» [1], where the limit of detection is calculated as  $(3 \cdot s_y)^*$ , with  $s_y$  as the residual standard deviation of a linear regression. With  $(10 \cdot s_y)^*$  the limit of quantification is calculated.

The linear working range was read out from a calibration curve.

The following table gives an overview over limit of detection and linear working range depending on the time for deposition.

Deposition	Limit of detection	Linear working range
30 s	2.1 µg/L	40 µg/L
60 s	0.9 µg/L	20 µg/L
90 s	0.7 µg/L	7 µg/L

### Comments

### Deposition

The deposition takes place when the method parameter «cleaning potential» is applied. When this potential is applied gas bubbles can be observed on the working electrode. The nascent hydrogen formed is necessary to reduce As(V) to As(III). The potential of -1 V is a compromise. On the one hand it has to be as negative as possible to generate enough nascent hydrogen for the reduction of As(V). On the other hand a too negative

potential results in worse reproducibility of the signals and reduced lifetime of the sensor.

### Adjustment of deposition potential

Different from the terminology in the software the deposition potential in this method is used for cleaning purposes. When this potential is applied other metals also deposited during the cleaning potential are stripped off the working electrode.

The deposition potential should be chosen as close as possible to the arsenic signal. Depending on the matrix the exact position of the arsenic signal can change. Therefore care has to be taken that the arsenic is not stripped off during this step.



Potential range for deposition potential

### Solution level



Make sure the electrode is fairly immersed into the measuring solution. A mark on the rear side of the electrode indicates the minimum (min) and maximum (max) solution level. If a measuring vessel other than 6.1415.150 is used volumes of sample and reagents may need to be adapted.

### Interference of organic substances

The influence of the organic substance on the peak of  $\beta(As)$  = 10 µg/L was tested.

EDTA	Does not interfere. Tested up to a concentration where EDTA starts to precipitate (approx. 0.03 mol/L).
Triton X-100	25 mg/L reduce the peak height of the As signal by approx 35%.

Organic interferences can be eliminated as described under «Sample preparation» for low and moderately contaminated samples.

### Interference of anions

The influence of the anion on the peak of  $\beta(As)$  = 10  $\mu g/L$  was tested.

I.	The whole curve shifts to higher currents. In the presence of c(iodide) = 0.1 mmol/L the peak height of the As signal is reduced by approx. 30%.
CI-	The chloride concentration has an influence on the reference potential. Therefore cleaning and deposition potential have to be adapted for samples with higher chloride concentration in order not to lose sensitivity. For sea water ( $c(Cl^-) \sim 0.5 \text{ mol/L}$ ), the potential shift is approx. +100 mV. Furthermore chloride facilitates the oxidation of gold. Therefore the positive potential in the conditioning cycles has to be adapted for higher chloride concentrations. For sea water, +0.85 V are used instead of +1.0 V.
NO <sub>2</sub> -	In the presence of 0.4 mol/L the peak height of the As signal is reduced by approx. 40%.

### Interference of cations

The influence of the cation on the peak of  $\beta(As) = 10 \ \mu g/L$  was tested. If nothing else is mentioned the interference was examined up to a concentration of 100 mg/L (equals 10000 time excess).

AI	Does not interfere.
Bi	Peak at +0.15 V, overlaps with the rear of the As signal. With 25 times excess the As signal vanishes in the Bi signal.
Са	Does not interfere.
Cd	Peak at -0.15 V, appears as shoulder on the front of the As signal. Up to 15 times excess do not show an influence on the peak height of the As signal.

Со	Does not interfere up to 15 times excess. With 1000 times excess the peak height of the As signal decreases by 50%.
Cr(VI)	Does not interfere up to 100 times excess. With 1000 times excess the peak height of the As signal decreases by 30% and an additional peak at +0.25 V appears.
Cu	Peak at +0.25 V, does not interfere up to about 5 times excess. With higher excess the Cu signal overlaps with the rear of the As signal and the sensitivity decreases. With 20 times excess the As signal vanishes in the Cu signal. Please see comment on «Reduction of copper interference».
Fe(III)	Does not interfere. With 100 times excess formation of a yellow Fe(III)- tri-citrate complex with the electrolyte and peak at +0.4 V.
Hg	Does not interfere up to 10 times excess. With 100 times excess formation of an amalgam working electrode which makes the determination of arsenic impossible.
Mg	Does not interfere.
Mn	Does not interfere up to 100 times excess. A 4000 times excess cause a deformation of the As signal and a reduced sensitivity. The As peak is approx. 30% smaller.
Ni	Does not interfere. Slightly reduced sensitivity for the As signal. With 4000 times excess peak at +0.25 V appears.
Pb	Peak at -0.15 V, does not interfere if equal in concentration. If higher in concentration it overlaps with the front of the As signal. With over 10 times excess also the sensitivity is affected. With 25 times excess the peak height of the As signal is reduced by approx. 60%.
Sb(III)	Peak at 0 V with the same sensitivity as arsenic.
Sb(V)	Does not interfere up to 1000 times excess. Higher concentrations show a peak at 0 V which appears as a

	shoulder on the As signal.
Se(IV)	Does not interfere if equal in concentration. But already a 5 times excess cause a deformation of the As signal and a loss in sensitivity. The As signal is about 30% smaller.
Se(VI)	Does not interfere up to 1000 times excess. Higher concentrations cause a deformation of the As signal and a loss in sensitivity.
Sn(IV)	10 times excess does not interfere. 100 times excess cause a deformation of the As signal and a peak at -0.2 V.
Zn	Does not interfere up to 10 times excess. Higher concentrations interfere with the formation of the nascent hydrogen which results in a reduced signal height. With 40 times excess the As signal is approx. 50%

### Reduction of copper interference

The concentration of copper and other cations in a sample can significantly be reduced using the IC-H sample preparation cartridge (6.1012.x10). Passing the sample through the IC-H cartridge cations are replaced by H<sup>+</sup> ions. In aqueous solution arsenic is present as arsenite ( $AsO_3^{3-}$ ) and arsenate ( $AsO_5^{3-}$ ), and is therefore not affected by the ion exchange.

### Procedure for the ion exchange:

For preparation the IC-H cartridge is rinsed with 10 mL ultrapure water. Then the sample is slowly forced through the cartridge using a disposable plastic syringe. The first few milliliters are discarded to avoid dilution of the sample. The cartridge can be regenerated by passing approx. 10 mL w(HCI) = 10% through.

### Storing the scTRACE Gold

When not in use, the electrode should be stored dry. After a series of determinations it is recommended to clean the electrode using the procedure described in the paragraph «Cleaning of the scTRACE Gold». Afterwards the electrode is rinsed thoroughly with ultrapure water.

Care has to be taken that the gold microwire is not damaged when the sensor is removed or stored outside the measuring stand.

# Method 2: Determination of arsenic(III)

With method 2 only the concentration of As(III) in the sample is determined. In this method the potential used for deposition is less negative than in method 1. At this potential only As(III) is reduced to As(0) and deposited on the gold working electrode. As(V) is not reduced and deposited under these conditions. During the subsequent stripping step the deposited As(0) is reoxidized to As(III) giving the analytical signal.

With 60 s for deposition the limit of detection is 0.3  $\mu$ g/L and the calibration is linear up to a concentration of 21  $\mu$ g/L.

### Sample preparation

Ground, drinking and mineral water can usually be analyzed directly.

If a digestion of the sample is necessary a quantification of As(III) is not possible. In this case only total arsenic can be determined, using method 1.

### Analysis

10 mL sample are pipetted into the measuring vessel and 2 mL electrolyte are added. The determination is carried out using the parameters given under «Parameters for As(III)».

The concentration of As(III) is quantified by two additions of As(III) standard solution.

### Measuring solution

10 mL.sample 2 mL electrolyte

### Parameters for As(III)

Determination	
No. of additions	2
No. of replications	2
Voltammetric	
Electrode	RDE/SSE
Measuring mode	SqW – Square wave
Stirring speed	2400 min <sup>-1</sup>
Hydrodynamic measurement	No
Conditioning cycles	
Start potential	-0.2 V
End potential	1 V
No. of cycles	10
Pretreatment	

Cleaning potential	-0.5 V
Cleaning time	60 s
Deposition potential	(-0.25 V) see «Comments»
Deposition time	5 s
Sweep	
Equilibration time	5 s
Start potential	-0.3 V
End potential	0.4 V
Potential step	0.01 V
Amplitude	0.02 V
Frequency	100 Hz
Sweep rate	1 V/s
Substance + calibration	
Calibration	Standard addition
Name	As(III)
Peak potential	0.0 V
Tolerance	0.05 V
Baseline	Linear
	Automatic

### Example for the determination of As(III)





### As(III) с = 1.436 ug/L +/-0.042 ug/L (2.90%) 500n 400n\_ 300n\_ <u></u> 200ŋ 100 -1.2e-006 -1.00u 1.00u 2.00u 3.00u

#### Results

Sample	β(As(III))
Bottled mineral water	1.4 µg/L*

c (q/L)

\* sample spiked

### Limit of detection and linear working range

The limit of detection was determined using the «regression approach» [1], where the limit of detection is calculated as « $3 \cdot s_y$ », with  $s_y$  as the residual standard deviation of a linear regression. With « $10 \cdot s_y$ » the limit of quantification is calculated.

The linear working range was readout from a calibration curve.

The following table gives an overview over limit of detection and linear working range depending on the time for deposition.

Deposition	Limit of	Linear working
	detection	range
30 s	0.5 µg/L	30 µg/L
60 s	0.3 µg/L	21 µg/L
90 s	0.2 µg/L	11 µg/L

### Comments

### Deposition

The deposition takes place when the method parameter «cleaning potential» is applied. This potential should not be more negative than -0.7 V otherwise also As(V) is partly reduced and deposited.

### Adjustment of deposition potential (potential 2)

See comments in method 1.

### Solution level

See comments in method 1.

### Interference of organic substances

See «Interference of organic substances» in method 1.

When As(III) should be determined it is not possible to eliminate organic interferences by digestion. A digestion would affect the oxidation state.

#### Interference of anions

-	Even traces of iodide make the
	determination impossible.

Other anions, see «Interference of anions» in method 1.

### Interference of cations

Interferences are only listed if they differ from method 1. Conditions as described in method 1.

Со	Does not interfere.
Cr(VI)	Oxidizes As(III) to As(V), therefore no As(III) can be expected when Cr(VI) is present.
Zn	Does not interfere.

### Reduction of copper interference

See comments in method 1.

Storing the scTRACE Gold

See comments in method 1.



# Cleaning of the scTRACE Gold

To maintain the scTRACE Gold the sensor can be electrochemically cleaned. The cleaning should be carried out before the sensor is stored and when it is put back into operation. But the cleaning can also be performed in between determinations to remove contaminations.

### Analysis

10 mL of the cleaning solution is pipetted into the measuring vessel. The cleaning is carried out using the parameters given under «Parameters for cleaning».

### Measuring solution

10 mL cleaning solution

### Parameters for cleaning

Determination	
No. of additions	0
No. of replications	4
Voltammetric	
Electrode	RDE/SSE
Measuring mode	DC – Sampled direct current
Stirring speed	2400 min <sup>-1</sup>
Hydrodynamic measurement	No
Conditioning cycles	
Start potential	-0.3 V
End potential	1 V
No. of cycles	10
Pretreatment	
Cleaning potential	0.1 V
Cleaning time	5 s
Deposition potential	-0.3 V
Deposition time	10 s
Sweep	
Equilibration time	5 s
Start potential	-0.3 V
End potential	0.2 V
Potential step	0.01 V
Potential step time	0.025 s
Sweep rate	0.4 V/s

### Example for the curves obtained during cleaning



### References

- [1] J. Mocak, A. Bond, S. Mitchell and G. Scollary, "A statistical overview of standard (IUPAC and ACS) and new procedures for determining the limits of detection and quantification: Application to voltammetric and stripping techniques," *Pure and Applied Chemistry*, vol. 69, no. 2, pp. 297-328, 1997.
- [2] P. Salaün, B. Planer-Friedrich and C. M. van den Berg, "Inorganic arsenic speciation in water and seawater by anodic stripping voltammetry with a gold microelectrode," *Analytica Chimica Acta*, vol. 585, no. 2, pp. 312-322, 2007.



## Appendix

### Method print for «Activation of the scTRACE Gold»

·	
Method parameters	
Method: AB416_ActivationTitle: Activation of thRemark1: 10 mL cleaning soRemark2: Cleaning solution	n scTRACE Gold.mth ne scTRACE Gold solution sn: 0.5 M H2SO4, 0.05 M KCl
Calibration : Standard additic Technique : Batch Addition : Manual	n
Sample ID : Activation scTRA Sample amount (mL): 10.000 Cell volume (mL): 10.000	ACE Gold
Voltammetric parameters	
Mode	: DC - Sampled Direct Current
Highest current range Lowest current range	: 10 mA : 100 nA
Electrode Stirrer speed (rpm)	: SSE/RDE : 2400
Initial electr. conditioning	: No
No. of additions No. of replications	: 0 : 4
Measure blank Addition purge time (s)	: No : 0
Initial purge time (s)	: 0
Conditioning cycles Start potential (V) End potential (V) No. of cycles	: -1.500 : 1.000 : 10
Hydrodynamic (measurement) Cleaning potential (V) Cleaning time (s) Deposition potential (V) Deposition time (s)	: No : 0.100 : 5.000 : -0.300 : 10.000
Sweep Equilibration time (s) Start potential (V) End potential (V) Voltage step (V) Voltage step time (s) Sweep rate (V/s)	$\begin{array}{cccccccccccccccccccccccccccccccccccc$
Cell off after measurement	: Yes
Peak evaluation	
Regression technique Peak evaluation Minimum peak width (V.steps) Minimum peak height (A) Reverse peaks Smooth factor Eliminate spikes	: Linear Regression : Height : 5 : 1.000e-010 : No : 4 : Yes
Substances Baseline	
Substance Addition automatic st	cart (V) end (V) type scope



### Method print «Method 1: Determination of total arsenic»

Method parameters							
Method       : AB416_Determination As(total) with scTRACE Gold.mth         Title       : Determination of As(total) with the scTRACE Gold         Remark1       : 10 mL sample + 2 mL electrolyte         Remark2       : Electrolyte: 1 M sulfamic acid, 0.5 M citric acid, 0.45 M KCl							
Calibration : Standar Technique : Batch Addition : Manual	ard addition 1						
Sample ID : sample Sample amount (mL): 10.00 Cell volume (mL): 12.000	0						
Voltammetric parameters							
Mode	: SqW - Square Wave						
Highest current range Lowest current range	: 10 mA : 100 nA						
Electrode Stirrer speed (rpm)	: SSE/RDE : 2400						
Initial electr. condition	ing : No						
No. of additions No. of replications	: 2 : 2						
Measure blank Addition purge time (s)	: No : 10						
Initial purge time (s)	: 10						
Conditioning cycles Start potential (V) End potential (V) No. of cycles	: -0.200 : 1.000 : 10						
Hydrodynamic (measurement Cleaning potential (V) Cleaning time (s) Deposition potential (V) Deposition time (s)	) : No : -1.000 : 60.000 : -0.250 : 5.000						
Sweep Equilibration time (s) Start potential (V) End potential (V) Voltage step (V) Amplitude (V) Frequency (Hz) Sweep rate (V/s)	: 5.000 : -0.300 : 0.400 : 0.010 : 0.020 : 100.000 : 0.992						
Cell off after measuremen	t : Yes						
Peak evaluation							
Regression technique Peak evaluation Minimum peak width (V.ste Minimum peak height (A) Reverse peaks Smooth factor Eliminate spikes	: Linear Regression : Height ps) : 5 : 1.000e-010 : No : 4 : Yes						
Substances							
As(total)	: 0.000 V +/- 0.050 V						
Standard solution Addition volume (mL)	: 2 1.000 mg/L : 0.050						
Total arsenic	: Final result (As(total)) = Conc * (12 / 10) * (1e+006 / 1) + 0 - 0						



### Method print «Method 2: Determination of arsenic(III)»

Method parameters							
Method: AB416_Determination As(III) with scTRACE Gold.mthTitle: Determination of As(III) with the scTRACE GoldRemark1: 10 mL sample + 2 mL electrolyteRemark2: Electrolyte: 1 M sulfamic acid, 0.5 M citric acid, 0.45 M KCl							
Calibration : Standar Technique : Batch Addition : Manual	libration : Standard addition chnique : Batch ddition : Manual						
Sample ID : As(III) Sample amount (mL): 10.00 Cell volume (mL): 12.000	) Mineral Water 30						
Voltammetric parameters							
Mode	: SqW - Square Wave						
Highest current range Lowest current range	: 10 mA : 100 nA						
Electrode Stirrer speed (rpm)	: SSE/RDE : 2000						
Initial electr. condition	ning : No						
No. of additions No. of replications	: 2 : 2						
Measure blank Addition purge time (s)	: No : 10						
Initial purge time (s)	: 10						
Conditioning cycles Start potential (V) End potential (V) No. of cycles	: -0.200 : 1.000 : 10						
Hydrodynamic (measurement Cleaning potential (V) Cleaning time (s) Deposition potential (V) Deposition time (s)	z) : No : -0.500 : 60.000 : -0.250 : 5.000						
Sweep Equilibration time (s) Start potential (V) End potential (V) Voltage step (V) Amplitude (V) Frequency (Hz) Sweep rate (V/s)	: 5.000 : -0.300 : 0.400 : 0.010 : 0.020 : 100.000 : 0.992						
Cell off after measuremen	nt : Yes						
Peak evaluation							
Regression technique Peak evaluation Minimum peak width (V.ste Minimum peak height (A) Reverse peaks Smooth factor Eliminate spikes	: Linear Regression : Height eps) : 5 : 1.000e-010 : No : 4 : Yes						
Substances	• 0.000 V +/- 0.050 V						
Standard solution Addition volume (mL)	: 2 1.000 mg/L : 0.020						
Arsenic(III)	: Final result (As(III)) = Conc * (12 / 10) * (1e+006 / 1) + 0 - 0						



### Method print «Cleaning of the scTRACE Gold»

Method parameters								
Method: AB416 Cleaning soTitle: Cleaning of the soRemark1: 10 mL cleaning soRemark2: Cleaning solution	TRACE Gold.mth scTRACE Gold plution 1: 0.5 M H2SO4, 0.05 M KCl							
Calibration : Standard addition Technique : Batch Addition : Manual	1							
Sample ID : Cleaning scTRACE Sample amount (mL): 10.000 Cell volume (mL): 10.000	Gold							
Voltammetric parameters								
 Mode	: DC - Sampled Direct Current							
Highest current range Lowest current range	: 10 mA : 100 nA							
Electrode Stirrer speed (rpm)	: SSE/RDE : 2400							
Initial electr. conditioning	: No							
No. of additions No. of replications	: 0 : 4							
Measure blank Addition purge time (s)	: No : 0							
Initial purge time (s)	: 0							
Conditioning cycles Start potential (V) End potential (V) No. of cycles	: -0.300 : 1.000 : 10							
Hydrodynamic (measurement) Cleaning potential (V) Cleaning time (s) Deposition potential (V) Deposition time (s)	: No : 0.100 : 5.000 : -0.300 : 10.000							
Sweep Equilibration time (s) Start potential (V) End potential (V) Voltage step (V) Voltage step time (s) Sweep rate (V/s)	: 5.000 : -0.300 : 0.200 : 0.010 : 0.025 : 0.397							
Cell off after measurement	: Yes							
Peak evaluation								
Regression technique Peak evaluation Minimum peak width (V.steps) Minimum peak height (A) Reverse peaks Smooth factor Eliminate spikes	: Linear Regression : Height : 5 : 1.000e-010 : No : 4 : Yes							
Substances								
Baseline								
Substance Addition automatic sta	art (V) end (V) type scope							

### Report «Example for the determination of As(total)»



### Report «Example for the determination of As(III)»



Sample ID Creator method Creator determ Modified by	: As(III) : p-fni : p-fni :	Mineral Water Date Date	e : 2013-05-0 e : 2013-09-1 e :	01 16	Time: Time: Time:	14:32:51 17:02:14	
Method Title Remark1 Remark2	: As(III) 5ppb.mth : As(III) Bestimmung : : Zitronensäure, Amidoschwefelsäure, KCl						
Sample amount Cell volume	: 10.000 : 12.000	mL mL					
Substance Conc. Conc.dev. Amount Add.amount	: As(III) : 1.196 : 0.035 : 14.356 : 20.000	ug/L ug/L (2.9 ng ng	90%)				
VR V	nA I.	mean Std.Dev.	I.delta	Comments			
1 - 1 0.017	143.0 1	42.0 1.356	5 0.0				
2 - 1 0.017	327.1 3	8.769	9 191.3				
$\begin{array}{cccccccccccccccccccccccccccccccccccc$	540.3 5 541.3	0.672	207.5				
Substance Ca	alibr.	Y.reg/offset	Slope	e Mean devi	at. Cor	r.Coeff.	
As(III) st	td.add.	1.416e-007	1.183e-003	1 6.102e-	009	0.99946	
Final results +/- Res. dev. % Comments							
As(III): Arsenic(III)	= 1.	436 ug/L (	0.042	2.896			