


# Complexometric titration of aluminum chloride using OMNIS system



With the OMNIS system a fast and accurate determination of aluminum in aluminum chloride by a complexometric back-titration with the copper ion-selective electrode is realized. The aluminum content is determined using cupric sulfate as titrant.

# Method description

## Sample

Aqueous solution of aluminum chloride

## Sample preparation

No sample preparation is required.

## Analysis

### Sample

1.4 mL sample solution is pipetted into the titration vessel and placed on the rack. Just before the titration, approximately 100 mL water, 5 mL acetate buffer and 10 mL EDTA solution are automatically added to the sample. After a reaction time of 3 min the solution is titrated with  $c(\text{CuSO}_4) = 0.1 \text{ mol/L}$  until after the equivalence point using the Cu-ISE.

## Configuration

|  |             |
|--|-------------|
| Main module Pick&Place S                                 | 2.1010.0010 |
| Pick&Place module  | 2.1014.0010 |
| "Peristaltic" (2-channel) pump module                    | 2.1016.0010 |
| Gripper fingers 42.8 - 65 mm                             | 6.02601.010 |
| Dummy panel for module plate                             | 6.02600.000 |
| OMNIS Rod Stirrer "Sample Robot"                         | 2.1006.0010 |
| Titration head 6xNS14 / 3xNS9 (P&P)                      | 6.01403.000 |
| Stirring propeller 30 mm ETFE                            | 6.01900.010 |
| OMNIS sample rack 9 x 250 mL, 2x                         | 6.02041.010 |
| Sample beaker (10x) PP 250 mL (P&P), 2x                  | 6.01400.100 |
| OMNIS Titrator (Advanced)                                | 2.1001.0210 |
| Cable MDL St/Bu 1 m, 3x                                  | 6.02102.020 |
| OMNIS Dosing Module, 2x                                  | 2.1003.0010 |
| OMNIS 5 mL cylinder unit, (acetate buffer)               | 6.03001.150 |
| OMNIS 10 mL cylinder unit, 2x (EDTA solution, titrant)   | 6.03001.210 |
| Analog measuring module                                  | 6.02101.010 |
| Cu-ISE   | 6.0502.140  |
| LL-ISE Reference electrode, Electrolyte c(KCl) = 3 mol/L | 6.0750.100  |
| Electrode cable plug-in head G / plug P, 1.5 m           | 6.02104.010 |
| Electrode cable, strand / 1 m / 2 x B                    | 6.2106.020  |

|   |             |
|---|-------------|
| OMNIS Stand-alone license (including one instrument license), OMNIS 1.0 | 6.06003.010 |
| OMNIS instrument license, 1x  | 6.06002.010 |

## Solutions

|                |   |
|----------------|---|
| Titrant        | $c(\text{CuSO}_4) = 0.1 \text{ mol/L}$ , if possible this solution should be bought from a supplier.          |
| Acetate buffer | Acetate buffer with $\text{pH} = 4.7$ , if possible this solution should be bought from a supplier.           |
| EDTA solution  | $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$ , if possible this solution should be bought from a supplier. |

## Parameters

|                   |  |
|-------------------|--|
| Mode              | MET U  |
| Pause             | 30 s   |
| Start volume      | $\text{Case}((10\text{-Sample Size}) > 0; 10\text{-Sample size}-1.0; 0)$ |
| Stirring rate     | 8  |
| Volume increment  | 50 $\mu\text{L}$   |
| Signal drift      | 30 mV/min  |
| Max. waiting time | 32 s   |
| Min. waiting time | 5 s  |
| Dosing rate       | Maximum  |
| Stop volume       | 10 mL  |
| Stop EP           | 1  |
| Volume after EP   | 1.0 mL   |
| EP criterion      | 5 mV   |
| EP recognition    | Greatest   |

## Results

|  |            |
|--|------------|
| Content $\beta_{\text{Al}^{3+}} / (\text{g/L})$ (n = 17) | s(rel) / % |
| 2.32   | 0.42       |

www.metrohm.com

