

Reliability for Routine pH Measurements



NEW

The proven Ecotrode now available as a maintenance-free gel electrode with aging indicator



Still within the permissible operating range: The decolorized zone of the gel gradually spreads up to the mark with use.

Highlights

- Ideal for **routine measurements** with similar types of samples
- **Maintenance-free gel electrolytes:** Refilling or checking the electrolyte is not necessary
- **Reliability through aging indicator:** The white gel loses its color with use.
As long as the decolorized zone is below the marking, accurate results are guaranteed.

Ordering information

6.0221.100	Ecotrode Gel (plug-in head G)
6.0221.600	Ecotrode Gel (plug-in head U)
6.2104.020	Electrode cable with plug-in head G / connector F
6.2104.600	Electrode cable with plug-in head U / connector F, 2 x B, 1 m

Visit us at www.metrohm.com