



Workshop

Sample preparation and particle characterisation made easy

You are invited

The Metrohm Australia team in Perth is inviting you to attend our showcase workshop.

Sample preparation and particle characterisation made easy.

Straight from the expert:

Matthew Johnson, our Retsch Product Manager for Australia and New Zealand, will guide you through the art of milling and particle sizing to also answer your questions first hand.

Gain valuable insights by watching Matthew perform live instrument demonstrations and teach you tips and tricks that only experts can teach.

Showcase equipment:

- The knife mill **GRINDOMIX GM 200** is ideal for grinding and homogenising foods and feeds. It allows to process sample volumes up to 0.7 litres quickly and reproducibly.
- The Ultra Centrifugal Mill **ZM 200** is a high speed rotor mill with unmatched performance. It is used for the rapid size reduction of soft to medium-hard and fibrous materials.
- The cyclone mill **TWISTER** is specially designed for the processing of foods and feeds for subsequent NIR analysis.
- The mixer mill **MM 400** is a compact versatile bench-top unit, which has been developed specially for dry, wet and cryogenic grinding of small amounts of sample.
- The vibratory sieve shakers of the series **AS 200** are used in quality control of raw materials, interim and finished products and in production monitoring.

Presenters:

- Matthew Johnson, Retsch Product Manager Australia and New Zealand
- Lavanya Avadiar, Retsch Territory Manager, Western Australia
- Penny Lewis, Senior Metrohm Territory Manager, Western Australia

15 May 2018

Western Australia

Metrohm Australia - Perth Office

Unit 6, 27 Tamara Drive
Cockburn Central WA 6164
Tel 08 6174 8200

Programme

09:30am - Arrival (Coffee)
09:45am - Presentation
10:30am - Morning Tea
10:45am - Live Demonstrations
10:30am - Q&A
12:00pm - End

Registration

Register online: mep.metrohm.com.au/mep-events/

For further information or enquiries feel free to contact **Paula Britz**: (+61) 02 8899 5222 or email paula.britz@metrohm.com.au

Metrohm Australia
Unit 11, 56 Buffalo Road, Gladesville, NSW 2111
PO Box 2173, Boronia Park, NSW 2111
T (+61) 02 8899 5200 F (+61) 02 8899 5299
mep.metrohm.com.au

Retsch[®]
MILLING SIEVING ASSISTING